

## Human CLDN18.2/PA-TU 8902 Stable Cell line

---

### I. Product Information:

<b>Catalog Number:</b>	C3001
<b>Cell Line Name:</b>	PA-TU 8902 stable cell clone expressing full-length human CLDN18.2 receptor
<b>Gene Synonyms:</b>	CLDN18.2
<b>Gene Sequence:</b>	Codon optimized from NP_001001026 (Met1-Val261)
<b>Protein Structure:</b>	Four span-transmembrane receptor
<b>Host Cell:</b>	PA-TU 8902, a human pancreatic adenocarcinoma cell line
<b>Quantity:</b>	Two vials of frozen cells (10 <sup>6</sup> cells/vial)
<b>Stability:</b>	>10 passages
<b>Freeze Media:</b>	90% FBS, 10% DMSO
<b>Storage:</b>	liquid nitrogen immediately upon receipt
<b>Culture Medium:</b>	DMEM with 10% FBS and 1ug/ml puromycin.
<b>Mycoplasma Test:</b>	Negative
<b>Application:</b>	Antibody binding assays, IHC/Western blot analysis, or use as cell immunogen

### II. Background:

Claudin-18 (CLDN18) is a member of a large family of four-span transmembrane proteins called Claudins. These proteins are the essential components of the mammalian tight junctions (TJs) in epithelial cells. Claudin-18 has two splice variants, 18.1 and 18.2. While CLDN18.1 is specifically expressed in the lung tissue, CLDN18.2 expression in normal tissue is more restricted and is only detected in small patches of stomach mucosal. CLDN18.2 expression is elevated in many types of epithelial cancers including stomach, esophagus, pancreatic and ovarian cancers. The expression of CLDN18.2 is not only detected in primary tumors, but also in the metastatic sites. Therefore, CLDN18.2 is an ideal target for monoclonal antibody-based cancer therapies.

---

*Disclaimer: For research use only. Not for use in humans.*




## Human CLDN18.2/PA-TU 8902 Stable Cell line

---

### III. Representative Data:

Detection of human CLDN18.2 expression on human CLDN18.2/PA-TU 8902 stable cells using a monoclonal antibody specific for human CLDN18.2 at two different concentrations (5 and 0.5ug/ml).



### IV. Thawing and Subculturing:

1. Remove the cell vial from liquid nitrogen tank and thaw cells quickly in a 37°C water bath.
2. In a T25 culture flask, add 10ml cell culture media which has been pre-warmed in a 37°C water bath for 30 minutes.
3. Transfer all the cells in the vial into the flask and place it inside of a CO<sub>2</sub> incubator with 5-8% CO<sub>2</sub>.
4. The cell doubling time should be 18-24 hours. Split the cells 1:5 or 1:10 every 3-5 days or when the cells are confluent.

### V. References:

Elsässer HP. *et al.* (1993): "Structural analysis of a new highly metastatic cell line Pa-Tu 8902 from a primary human pancreatic adenocarcinoma". *Virchows Arch B Cell Pathol Incl Mol Pathol.* 64(4):201-7.

Türeci O. *et al.* (2011): "Claudin-18 gene structure, regulation, and expression is evolutionary conserved in mammals". *Gene*, 481(2), p83-92.

Sahin U. *et al.* (2008): "Claudin-18 Splice Variant 2 Is a Pan-Cancer Target Suitable for Therapeutic Antibody Development". *Clin. Cancer Res.* 14 (23) p7624-7634.

---

*Disclaimer: For research use only. Not for use in humans.*



## Human CLDN18.2/PA-TU 8902 Stable Cell line

---

Niimi T. *et al.* (2001): "claudin-18, a Novel Downstream Target Gene for the T/EBP/NKX2.1 Homeodomain Transcription Factor, Encodes Lung- and Stomach-Specific Isoforms through Alternative Splicing". *Mol. Cell. Biol.* 21(21), p7380-7390.

### VI. Limited Usage License Agreement:

This is a legal agreement between the buyer (Licensee) and Accurus Biosciences Inc. regarding the use of Accurus Biosciences' stable cell line products. By purchasing and using the stable cell lines, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- A) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owner by the buyer or to a subcontractor will be permitted only when written approval from Accurus Biosciences is obtained.
- B) The buyer cannot sell or otherwise transfer a) this product, b) its components or c) materials made using this products to a third party.
- C) The products sold by Accurus Biosciences are for laboratory and animal research use only. The products are not to be used on human studies or consumption.

---

*Disclaimer: For research use only. Not for use in humans.*